



## IgG Product Specifications:

IgG is the first test we released on our Sidekick platform. Sometimes you can't wait for the lab. When a foal's life depends on timely results, Sidekick gives you real-time data for an accurate diagnosis and on-the-spot treatment.

Use Sidekick IgG tests as a routine protocol for screening foals that present an elevated risk of failure of passive transfer.

An accurate measurement of IgG levels is essential to ensure treatment is initiated when passive transfer is unsuccessful.

Sidekick gives easy to read and fully quantitative results between 80-2000mg/dL within 5 minutes, which allows you to test twice or more within 24 hours to monitor progress.

